



Specification

8.00028.0100 Acetophenone for synthesis

Specification		
Colour (visual)	colorless to yellowish	
Appearance of substance (visual)	clear liquid	
Assay (GC, area%)	≥ 99.0	% (a/a)
Density (d 20 °C/ 4 °C)	1.025 - 1.028	
Identity (IR)	passes test	

Due to its specific melting range the product may be solid, liquid, a solidified melt or a supercooled melt.

Dr. Oliver Schramel
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.